



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® fan heater with 2 power/ 3 heat settings and cold switch, 2100 W

**€42.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Handy heater fan with various heat and speed combinations for everything that needs to be dried in the laboratory.  
Max. Heating approx. 2100 W, overheating protection.  
2 power/3 heat levels, additional cold air button.  
Housing with soft-touch surface, suspension eyelet.  
**Color:** silver

## Additional Information

---

No.	1-7015
Manufacturer (Brand)	neoLab
EAN	4058072086565
Color	gray
Voltage	230 V
Heaters and warming cabinets features	with thermal protection shutdown
Number of blower stages	2
Number of heating stages	3
Functionality	electric
Power	2100 W

