



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Corning Shaker Incubator, 49L, 230V EU&UK Power Plug

€5,089.28
plus VAT &
Shipping

Product Images



Description

Stackable incubation shaker with door at the front and compact footprint, for maximum 4 x 1L culture flasks. A wide range of platforms and clamps are available for this incubation shaker for mounting flasks, Petri dishes and test tubes.

Uniform temperature distribution, wide temperature range, precise control of shaking motion.

Volume: 49 liters

Additional Information

No.	CQ-0198
Manufacturer (Brand)	Corning
old neoLab article no.	C-2061
EAN	4058072225759
Transport temperature	Room temperature
Color	White Grey
Temp MAX	80 °C
Temp MIN	20 °C
Volume	49 l
Wide	550 mm
Width/inside	375 mm
Height	580 mm
Height inside	381.0 mm
Length	425 mm
Depth	425 mm
Depth inside	343 mm
Number of shelves/ grids	1
Operation via	Membrane keypad
Weight	37.7 kg
Length inside	343 mm
Voltage	230 V
Heaters and warming cabinets features	with overheating protection with infinitely variable temperature control with integrated shaker with timer with glass door
Power	780 W
Type of convection	Natural
Temp. accuracy	+/- 0.1 °C
Type incubator/ incubator	without CO2 supply
Material inside	Stainless steel

Page door suspension	right
Positioning variants	single
Shaker type	Orbital/circular shaker
Type temperature control	PID
Type tempering cabinet	Incubator with integrated shaker

