



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Corning Shaking Incubator, 71I, 230V EU&UK power plug

**€5,796.15  
plus VAT &  
Shipping**

### Product Images



## Description

---

Incubation shaker with door at the front and compact footprint, for a maximum of 4 x 2 liter culture flasks. A wide range of platforms and clamps are available for this incubation shaker for mounting flasks, Petri dishes and test tubes.

Uniform temperature distribution, large temperature range, precise control of shaking motion.

**Volume: 71 liters**

## Additional Information

No.	CQ-0200
Manufacturer (Brand)	Corning
old neoLab article no.	C-2068
EAN	4058072225766
Transport temperature	Room temperature
Color	White Grey
Temp MAX	80 °C
Temp MIN	20 °C
Volume	71 l
Wide	544 mm
Width/inside	396 mm
Height	635 mm
Height inside	406.0 mm
Length	574 mm
Depth	574 mm
Depth inside	442 mm
Number of shelves/ grids	1
Operation via	Membrane keypad
Weight	63.6 kg
Length inside	442 mm
Voltage	230 V
Heaters and warming cabinets features	with overheating protection with infinitely variable temperature control with integrated shaker with timer with glass door
Power	920 W
Type of convection	Natural
Temp. accuracy	+/- 0.1 °C
Type incubator/ incubator	without CO2 supply
Material inside	Stainless steel

Page door suspension	right
Positioning variants	single
Shaker type	Orbital/circular shaker
Type temperature control	PID
Type tempering cabinet	Incubator with integrated shaker

