



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Zefa STACKABLE SHELL INCUBATORS – PREMIUM, SOK2190, 2x190 liters

**Price on  
request  
plus VAT &  
Shipping**

### Product Images

## Description

---

The SOK series **shaker incubators** are stackable. The capacity is tripled. This means that you only need the space for one incubator for three stacked units. All models feature an insulated, fold-down door with a double-layered glass window so you can easily view the interior. On all cooling models, microprocessor control provides unmatched versatility for the user to create a customized program (up to 9 segments, with cycles) to automate changes in functional parameters.

- Stackable up to three units for maximum space savings
- 6" LCD touch screen panel
- Direct beltless drive system is durable, economical, maintenance-free and allows an extended speed range
- from 30 to 300 rpm, + 1 rpm with minimal vibration
- Wide-angle hinged door with door handle and an ergonomically extendable shaking platform provide convenient access to your samples
- The top of the incubator can be further used as a work area or as a shelf for small labware

## Additional Information

No.	8-8202
Manufacturer (Brand)	Zefa
EAN	4066292447653
Transport temperature	Room temperature
Temp MAX	60 °C
Temp MIN	4 °C
Volume	190000 ml
Width/inside	920 mm
Height inside	395.0 mm
Depth inside	532 mm
Number of shelves/ grids	2
Display (reading, setting)	digital
Weight	430 kg
Voltage	240 V
Heaters and warming cabinets features	with infinitely variable temperature control with integrated shaker with timer with alarm function programmable
Power	2400 W
Loading	from the front
Temp. accuracy	+/- 0.1 °C
Type incubator/ incubator	without CO2 supply
for vessels with volume up to	2000 ml
Shaker type	Orbital/circular shaker

