



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Binder Series B - Incubators with mechanical adjustment, 28 L

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Advantages

- Safe and reproducible results even with high batch throughputs in long-term operation
- Inner glass door for a stable atmosphere

Important features

- Temperature range: 30 °C to 70 °C
- Hydraulic-mechanical thermostat
- Inner door made of safety glass (ESG)
- 2 chrome-plated slide-in grids
- stackable

Additional Information

No.	BB-1185
Manufacturer (Brand)	Binder
EAN	4058072522001
Transport temperature	Room temperature
Color	gray red
Temp MAX	70 °C
Temp MIN	30 °C
Wide	580 mm
Width/inside	400 mm
Height	405 mm
Height inside	280.0 mm
Depth	425 mm
Depth inside	250 mm
Number of shelves/ grids	2
Display (reading, setting)	analog
Operation via	Knobs/Buttons
Weight	23 kg
Voltage	230 V
Heaters and warming cabinets features	with interior glass door
Functionality	mechanical
Power	250 W
Type of convection	Natural
Loading	from the front
Type incubator/ incubator	without CO2 supply
Type tempering cabinet	Incubator/ incubator

