



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Binder Series BD - Incubators Avantgarde.Line with free convection, 260 L

€3,763.00
plus VAT &
Shipping

Product Images



Description

The BD Avantgarde.Line series is perfectly suited for gentle incubation of organisms and for conditioning heat-sensitive media. The new unit design and the homogeneous temperature distribution ensure best results.

Advantages

- Very good temporal and spatial temperature accuracy
- High process reliability
- 100 °C disinfection routine

Important features

- Temperature range: +5 °C above ambient to +100 °C
- Up to 30% more energy efficient than commercially available devices
- High temperature accuracy due to APT.line™ technology
- Natural convection
- Controller with LCD display
- Electromechanical control of the exhaust air flap
- Inner door made of safety glass (ESG)
- 2 chrome-plated slide-in grids
- Stackable units up to 115 L
- Integrated temperature selector monitor class 3.1 (DIN 12880) with visual alarm
- Ergonomic handle design
- USB port for data recording

Additional Information

No.	BB-1247
Manufacturer (Brand)	Binder
EAN	4058072522629
Transport temperature	Room temperature
Color	red grey
Temp MAX	100 °C
Wide	810 mm
Width/inside	610 mm
Height	880 mm
Height inside	760.0 mm
Depth	760 mm
Depth inside	545 mm
Number of shelves/ grids	2
Display (reading, setting)	digital
Weight	83 kg
Voltage	230 V
Heaters and warming cabinets features	Temp MIN/MAX depends on ambient temperature with internal glass door
Type of convection	Natural
Loading	from the front
Power	850 W
Model incubator/ incubator	BINDER BD 240/260
Temp MIN above ambient temperature	5 °C
Temp. accuracy	+/- 0.3 °C
Type incubator/ incubator	without CO2 supply
Type tempering cabinet	Incubator/ incubator

