



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Binder Series BD - Incubators Classic.Line with free convection, 400 L

**€5,383.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The BINDER BD series incubator with natural convection is a specialist for long-term stable continuous use. It is perfectly suited for gentle incubation of organisms and for conditioning heat-sensitive media.

### Advantages

- Reliable and reproducible results even with high batch throughputs in long-term operation
- Disinfection routine at 100 °C
- Internal glass door for a stable atmosphere

### Important features

- Temperature range: +5 °C above ambient to +100 °C
- APT.line™ preheating chamber technology
- Natural convection
- Adjustable exhaust air flap
- Controller with time functions
- Inner door made of safety glass (ESG)
- 2 chrome-plated slide-in grids
- Stackable units up to 115 L
- Temperature selector monitor class 3.1 (DIN 12880) with optical temperature alarm
- Computer interface: RS 422

## Additional Information

No.	BB-1193
Manufacturer (Brand)	Binder
EAN	4058072522087
Transport temperature	Room temperature
Color	red grey
Temp MAX	100 °C
Wide	1235 mm
Width/inside	1000 mm
Height	1025 mm
Height inside	800.0 mm
Depth	765 mm
Depth inside	500 mm
Number of shelves/ grids	2
Display (reading, setting)	digital
Weight	135 kg
Voltage	120 V
Heaters and warming cabinets features	with internal glass door Temp MIN/MAX depends on ambient temperature
Power	850 W
Type of convection	Natural
Loading	from the front
Temp. accuracy	+/- 0.5 °C
Type incubator/ incubator	without CO2 supply
Model incubator/ incubator	BINDER BD 400
Temp MIN above ambient temperature	5 °C
Type tempering cabinet	Incubator/ incubator

