



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Binder Series BD-S Solid.Line - Incubators with free convection, 62 L

€1,294.00
plus VAT &
Shipping

Product Images



Description

Series BD-S Solid.Line | Standard Incubators with Natural Convection

The new products of the Solid.Line convince with proven BINDER quality and reliability. They are suitable for incubation. The performance features of the Solid.Line are geared to the respective application and are thus used in research and quality assurance.

Advantages

- Easy to operate
- Precise temperature control
- Easy and ergonomic door opening

Important features

- Temperature range: +5 °C above ambient to +70 °C
- APT.line™ preheating chamber technology
- Natural convection
- Adjustable exhaust air flap
- Controller with timer function
- Inner door made of safety glass (ESG)
- 1 chrome-plated slide-in grille, incl. shelf support
- Integrated temperature selection monitor class 3.1 (DIN 12880) with optical alarm

Additional Information

No.	BB-1657
Manufacturer (Brand)	Binder
EAN	4058072524081
Transport temperature	Room temperature
Color	gray red
Temp MAX	70 °C
Wide	614 mm
Width/inside	400 mm
Height	625 mm
Height inside	440.0 mm
Depth	577 mm
Depth inside	350 mm
Heaters and warming cabinets features	Temp MIN/MAX depends on ambient temperature with internal glass door with stepless temperature control
Number of shelves/ grids	1
Display (reading, setting)	digital
Type of convection	Natural
Operation via	Knobs/Buttons
Loading	from the front
Weight	36 kg
Model incubator/ incubator	BINDER BD 53/56
Useful volume	62000 ml
Voltage	230 V
Temp MIN above ambient temperature	5 °C
Type incubator/ incubator	without CO2 supply

