



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Binder Series BD-S Solid.Line - Incubators with free convection, 118 L

**€1,757.00**  
**plus VAT &**  
**Shipping**

### Product Images



# Description

---

## Series BD-S Solid.Line | Standard Incubators with Natural Convection

The new products of the Solid.Line convince with proven BINDER quality and reliability. They are suitable for incubation. The performance features of the Solid.Line are geared to the respective application and are thus used in research and quality assurance.

### Advantages

- Easy to operate
- Precise temperature control
- Easy and ergonomic door opening

### Important features

- Temperature range: +5 °C above ambient to +70 °C
- APT.line™ preheating chamber technology
- Natural convection
- Adjustable exhaust air flap
- Controller with timer function
- Inner door made of safety glass (ESG)
- 1 chrome-plated slide-in grille, incl. shelf support
- Integrated temperature selection monitor class 3.1 (DIN 12880) with optical alarm

## Additional Information

No.	BB-1664
Manufacturer (Brand)	Binder
EAN	4058072524159
Transport temperature	Room temperature
Color	gray red
Temp MAX	70 °C
Wide	764 mm
Width/inside	550 mm
Height	735 mm
Height inside	550.0 mm
Depth	616 mm
Depth inside	390 mm
Number of shelves/ grids	1
Display (reading, setting)	digital
Operation via	Knobs/Buttons
Weight	47 kg
Voltage	120 V
Heaters and warming cabinets features	Temp MIN/MAX depends on ambient temperature with internal glass door with stepless temperature control
Type of convection	Natural
Useful volume	118000 ml
Loading	from the front
Type incubator/ incubator	without CO2 supply
Model incubator/ incubator	BINDER BD 115
Temp MIN above ambient temperature	5 °C

