



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Binder Series BF – Incubators Avantgarde.Line with circulating air, 260 L

**€4,557.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### **BF Avantgarde.Line series | Standard incubators with circulating air**

The BINDER BF series incubator with convection is suitable for all incubation tasks, especially when fully loaded and for large batch throughputs. The latest APT.line technology provides absolutely homogeneous temperature distribution. The extended time functions and adjustable fan allow ideal control of temperature parameters and recirculation conditions.

#### Advantages

- Best energy efficiency on the market
- Very good temporal and spatial temperature accuracy
- High process reliability
- 100 °C disinfection routine

#### Important features

- Temperature range: +7–8 °C above ambient to +100 °C
- Up to 30% more energy efficient than commercially available devices
- High temperature accuracy due to APT.line™ technology
- Circulating air
- Controller with LCD display
- Electromechanical control of the exhaust air flap
- Inner door made of safety glass (ESG)
- 2 chrome-plated slide-in grids
- Stackable units up to 115 L
- 4 sturdy castors, two with locking brake, from 720 L
- Integrated temperature selector monitor class 3.1 (DIN 12880) with visual alarm
- Ergonomic handle design
- USB port for data recording

## Additional Information

No.	BB-1239
Manufacturer (Brand)	Binder
EAN	4058072522544
Transport temperature	Room temperature
Color	gray red
Temp MAX	100 °C
Wide	810 mm
Width/inside	650 mm
Height	845 mm
Height inside	780.0 mm
Depth	760 mm
Depth inside	510 mm
Number of shelves/ grids	2
Display (reading, setting)	digital
Operation via	Knobs/Buttons
Weight	82 kg
Voltage	230 V
Heaters and warming cabinets features	with inside glass door Temp MIN/MAX depends on ambient temperature with infinitely variable temperature control
Type of convection	Circulating air
Loading	from the front
Model incubator/ incubator	BINDER BF 240/260
Useful volume	257000 ml
Temp MIN above ambient temperature	7 °C
Temp. accuracy	+/- 0.2 °C
Type incubator/ incubator	without CO2 supply

