



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Binder Series BF – Incubators Classic.Line with circulating air, 400 L

€5,696.00
plus VAT &
Shipping

Product Images



Description

BF Classic.Line series | Standard incubators with circulating air

The BINDER BF series incubator with circulating air is suitable for all gentle incubation tasks, especially when fully loaded and for large batch throughputs. This incubator operates extremely homogeneously and has fast recovery times.

Advantages

- Reliable and reproducible results even with high batch throughputs in long-term operation
- Disinfection routine at +100 °C
- Inner glass door for a stable atmosphere

Important features

- Temperature range: +5 °C above ambient to +100 °C
- APT.line™ preheating chamber technology
- Adjustable fan speed
- Adjustable exhaust air damper
- Controller with time functions
- Inner door made of safety glass (ESG)
- 2 chrome-plated slide-in grids
- Temperature selector monitor class 3.1 (DIN 12880) with optical temperature alarm
- Computer interface: RS 422

Additional Information

No.	BB-1209
Manufacturer (Brand)	Binder
EAN	4058072522247
Transport temperature	Room temperature
Color	gray red
Temp MAX	100 °C
Wide	1235 mm
Width/inside	1000 mm
Height	1025 mm
Height inside	800.0 mm
Depth	765 mm
Depth inside	500 mm
Number of shelves/ grids	2
Display (reading, setting)	digital
Operation via	Knobs/Buttons
Weight	145 kg
Voltage	230 V
Heaters and warming cabinets features	with infinitely variable temperature control with interior glass door Temp MIN/MAX is dependent on ambient temperature
Type of convection	Circulating air
Useful volume	400000 ml
Loading	from the front
Temp. accuracy	+/- 0.4 °C
Type incubator/ incubator	without CO2 supply
Model incubator/ incubator	BINDER BF 400
Temp MIN above ambient temperature	5 °C

