



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wiggens Shaking Incubator, 53 liters, 15 °C to 60 °C

€5,024.92
plus VAT &
Shipping

Product Images



Description

Combining a benchtop incubator and a double-acting shaker (orbital and reciprocating), **Wiggens shaker incubators** offer versatility for a variety of biological and biotechnological applications. Ideal for bacterial suspensions, hybridizations, staining and destaining, metabolic studies, and for culturing cells

- Dual shaking modes: recirculation or reciprocation
- Speed range: 10 to 300 rpm
- Programmable shaking motion
- Timer range: 10 sec. to 999 hr. 59 min. 59 sec
- Forward & reverse, pause: 10 sec. to 59 min. 59 sec
- Selectable orbit diameter or stroke length: 20, 30, 40 mm

Delivery without attachment/platform !

Additional Information

No.	WI-1560
Manufacturer (Brand)	Wiggins
old neoLab article no.	F-1043
EAN	4058072163648
Transport temperature	Room temperature
Color	Blue. White
Material	Stainless steel Painted steel
Temp MAX	60 °C
Temp MIN	25 °C
Volume	53 l
Wide	440 mm
Width/inside	410 mm
Height	625 mm
Height inside	320.0 mm
Depth	740 mm
Depth inside	410 mm
Display (reading, setting)	digital
Operation via	Rotary knobs/Keys Foil keyboard
Weight	85 kg
Voltage	230 V
Type interface	RS232
Heaters and warming cabinets features	with timer
Resolution °C	0.1 °C
Functionality	electronic
Heating capacity	600 W
Type of convection	Natural
Loading	from the front

Temp. accuracy	+/- 0.1 °C
Type incubator/ incubator	without CO2 supply
Cooling with	Air
Model incubator/ incubator	Wiggens WS-300R
Positioning variants	single
Shaker type	Orbital/circular shaker
Type temperature control	PID
Type tempering cabinet	Incubator with integrated shaker

