



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Muffle Furnace, 6.7 L, 1300 °C

€3,723.00
plus VAT &
Shipping

Product Images



Description

The **neolab muffle furnaces** are suitable for thermal processes for materials up to 1300 °C: Drying, heating, curing, cremating. The muffle furnaces can be used in scientific laboratories, schools and other educational institutes, in medicine and in industry. The chambers are made of thermally robust ceramic with fiberglass insulation.

- heat-resistant ceramic materials
- short heating times
- Thermoregulator
- high degree of accuracy

Additional functions can be added as an option

- Process monitoring window up to 1100 °C
- Chimney with fan
- Replacement ceramic shelves
- Device for protective gas atmosphere in the furnace
- programmable temperature controller
- Digital timer
- Buzzer

Additional Information

| | |
|----------------------------|------------------|
| No. | 7-9352 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072195229 |
| Transport temperature | Room temperature |
| Color | light gray |
| Material | Steel |
| Temp MAX | 1300 °C |
| Temp MIN | 10 °C |
| Volume | 6.7 l |
| Wide | 445 mm |
| Width/inside | 145 mm |
| Height | 525 mm |
| Height inside | 135.0 mm |
| Length | 575 mm |
| Depth | 575 mm |
| Depth inside | 310 mm |
| Display (reading, setting) | digital |
| Operation via | Knobs/Buttons |
| Weight | 35 kg |
| Length inside | 310 mm |
| Voltage | 230 V |
| Number of power levels | 1 |
| Power | 2400 W |
| Temp. accuracy | +/- 1.0 °C |
| Type tempering cabinet | Muffle furnace |

