



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoBlock Heater Duo, without heating block

Special Price  
**€929.00** was  
**€1,100.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Thermoblock heater for 2 aluminum blocks in a row or for one large block for microtiter plates.  
Small footprint, great technical advantages:

- the stored values are retained even in the event of a power failure.
- the temperature probe can be inserted from the top or from the back of the unit.
- easy temperature readjustment.
- 11 mm high digital display for actual and set temperature.
- Residual heat indication by flashing light when temperature exceeds 60°C.
- Timer for automatic switch-off.
- additional power switch on the back.

### **Without heating block**

#### Technical data:

Timer: 0-99.9 h

Residual heat indicator above 60 °C

Temperature: RT to 110°C

Resolution: 0.1°C

Connection: 230 Volt / 50 Hz

Power consumption: max. 160 VA

Dimensions: 137 x 184 x 151 mm

Weight: 2.5 kg

## Additional Information

---

No.	2-2504
Manufacturer (Brand)	neoLabLine
EAN	4058072061142
Transport temperature	Room temperature
Temp MAX	110 °C
Wide	140 mm
Height	150 mm
Length	190 mm
Depth	190 mm
Operation via	Membrane keypad
Weight	2.6 kg
Voltage	230 V
Resolution °C	0.1 °C
Energy supply	Power grid
Power	200 W
Temp MIN heater	Room temperature

