



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLabLine neoBlock Heater Mono I, without block

Special Price  
**€649.00** was  
**€799.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Low-cost heating block thermostat with great technical advantages and small footprint.

The data programmed in the microprocessor remain stored even in case of power failure.

The temperature sensor can be inserted from the top or from the back of the unit.

Large, 11 mm high digital display for temperature. Simultaneous display of set and actual temperature. The unit can be readjusted at any time. When the block is hotter than 60°C and the unit is switched off, a flashing light warns of residual heat.

Timer for automatic heater shutdown adjustable from 0-99.9 hours.

An additional power switch on the back provides more safety.

**Without heating block.**

Technical data:

Timer: 0-99.9 h

Residual heat indicator at temperatures above 60°C

Temperature: RT to 110°C

Resolution : 0.1°C

Connection: 230 volts /50 Hz

Power consumption: max. 160 watts

Weight: 1.8 kg

Dimensions : 107 x 143 x 151 mm

## Additional Information

---

No.	2-2503
Manufacturer (Brand)	neoLabLine
EAN	4058072061135
Transport temperature	Room temperature
Temp MAX	110 °C
Wide	107 mm
Height	151 mm
Length	143 mm
Depth	143 mm
Operation via	Membrane keypad
Weight	1.8 kg
Voltage	230 V
Resolution °C	0.1 °C
Energy supply	Power grid
Power	160 W
Temp MIN heater	Room temperature

