



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Corning Table Shaker Incubator W/Platform, 230V EU&UK Power Plug

**€3,675.61**  
plus VAT &  
Shipping

### Product Images



## Description

---

Compact benchtop incubation shaker for a maximum of 4 x 1 liter culture flasks. The included software records speed and temperature. Uniform temperature distribution, wide temperature range, precise control of shaking motion.

## Additional Information

No.	CQ-0234
Manufacturer (Brand)	Corning
old neoLab article no.	C-2051
EAN	4058072225780
Transport temperature	Room temperature
Color	White Grey
Temp MAX	70 °C
Temp MIN	20 °C
Volume	10 l
Wide	370 mm
Width/inside	360 mm
Height	400 mm
Height inside	300.0 mm
Length	530 mm
Depth	530 mm
Depth inside	520 mm
Number of shelves/ grids	1
Operation via	Membrane keypad
Weight	20.2 kg
Length inside	520 mm
Voltage	230 V
Type interface	USB
Heaters and warming cabinets features	with overheating protection with infinitely variable temperature control with integrated shaker with timer with interface
Power	700 W
Type of convection	Natural
Temp. accuracy	+/- 0.1 °C
Type incubator/ incubator	without CO2 supply

Material inside	Polycarbonate (PC)
Positioning variants	single
Shaker type	Orbital/circular shaker
Type temperature control	PID
Type tempering cabinet	Incubator with integrated shaker

