



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Kettle plastic housing, 1.5 l, 2200 W

€43.70
plus VAT &
Shipping

Product Images



Description

Boiling water within a very short time.

Both stoves are cordless through separate device base and have:

- concealed heating element
- Over-temperature protection
- Automatic steam stop
- Pilot light
- Cable rewind

Connection: 230 Volt / 50 Hz

Order No. 2-1162

Capacity: 1.5 l

Heating capacity: 2200 W

transparent water level indicator on both sides

retention sieve

Additional Information

No.	2-1162
Manufacturer (Brand)	neoLab
EAN	4058072055707
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	High-density polyethylene (HDPE)
Volume	1500 ml
sterile	No
Height	220 mm
Operation via	Knobs/Buttons
Weight	0.9 kg
Voltage	230 V
Power	2200 W

