



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wessamat flake ice maker with air cooling, capacity 40 kg/day

**€3,300.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

Compact, versatile ice maker (ice machine) with air cooling for the production of flake ice (temperature  $-0.5^{\circ}\text{C}$ ). Ideally suited for use in laboratories and hospitals. Tailored to users who need a space-saving and, above all, cost-effective device for the production of flake ice. Stocking takes place as standard in the integrated storage container.

All units are equipped with the environmentally friendly refrigerant R290 (propane).

## Additional Information

No.	1-0001
Manufacturer (Brand)	Wessamat
EAN	4058072115586
Color	silver
Material	Stainless steel
Temp MAX	-0.5 °C
Temp MIN	-0.5 °C
Wide	520 mm
Height	637 mm
Length	345 mm
Depth	345 mm
Operation via	Knobs/Buttons
Weight	39 kg
Refrigerant medium	Refrigerant (closed system)
Refrigerant closed system	R 290 (propane)
Voltage	230 V
for ambient temperature	10 to 30°C
Refrigerators and freezers features	suitable for use with deionized water
Ice capacity kg/24h	40 kg/24h
Quantity compressors	1
Capacity reservoir	10 kg
Power	410 W
Devices self cooling	Air cooling
Positioning variants	single
Type refrigerator / freezer	Flake ice maker

