



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wessamat flake ice maker with air cooling, capacity 80 kg/day

€4,090.00
plus VAT &
Shipping



Description

Compact, versatile ice maker (ice machine) with air cooling for the production of flake ice (temperature -0.5°C). Ideally suited for use in laboratories and hospitals. Tailored to users who need a space-saving and, above all, cost-effective device for the production of flake ice. Stocking takes place as standard in the integrated storage container.

All units are equipped with the environmentally friendly refrigerant R290 (propane).

Additional Information

No.	1-0002
Manufacturer (Brand)	Wessamat
EAN	4058072115609
Color	silver
Material	Stainless steel
Temp MAX	-0.5 °C
Temp MIN	-0.5 °C
Wide	586 mm
Height	850 mm
Length	485 mm
Depth	485 mm
Operation via	Knobs/Buttons
Weight	55 kg
Refrigerant medium	Refrigerant (closed system)
Refrigerant closed system	R 290 (propane)
Voltage	230 V
for ambient temperature	10 to 30°C
Refrigerators and freezers features	suitable for use with deionized water
Ice capacity kg/24h	80 kg/24h
Quantity compressors	1
Capacity reservoir	25 kg
Power	520 W
Devices self cooling	Air cooling
Positioning variants	single
Type refrigerator / freezer	Flake ice maker

