



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Wessamat flake ice maker with air cooling,
capacity 135 kg/day**

**€4,726.00
plus VAT &
Shipping**

Product Images



Description

Compact, versatile ice maker (ice machine) with air cooling for the production of flake ice (temperature -0.5°C). Ideally suited for use in laboratories and hospitals. Tailored to users who need a space-saving and, above all, cost-effective device for the production of flake ice. Stocking takes place as standard in the integrated storage container.

On request: units also available with environmentally friendly refrigerant R290 (propane).

Additional Information

No.	1-0003
Manufacturer (Brand)	Wessamat
EAN	4058072115623
Color	silver
Material	Stainless steel
Temp MAX	-0.5 °C
Temp MIN	-0.5 °C
Wide	588 mm
Height	1020 mm
Length	680 mm
Depth	680 mm
Operation via	Knobs/Buttons
Weight	69 kg
Refrigerant medium	Refrigerant (closed system)
Refrigerant closed system	R 290 (propane)
Voltage	230 V
for ambient temperature	10 to 30°C
Refrigerators and freezers features	suitable for use with deionized water
Ice capacity kg/24h	135 kg/24h
Quantity compressors	1
Capacity reservoir	42 kg
Power	590 W
Devices self cooling	Air cooling
Positioning variants	single
Type refrigerator / freezer	Flake ice maker

