



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wessamat flake ice maker with water cooling, capacity 40 kg/day

€3,385.00
plus VAT &
Shipping

Product Images



Description

Compact, versatile ice maker (ice machine) with air cooling for the production of flake ice (temperature -0.5°C). Ideally suited for use in laboratories and hospitals. Tailored to users who need a space-saving and, above all, cost-effective device for the production of flake ice. Stocking takes place as standard in the integrated storage container.

All units are equipped with the environmentally friendly refrigerant R290 (propane).

Additional Information

No.	1-0004
Manufacturer (Brand)	Wessamat
EAN	4058072115593
Color	silver
Material	Stainless steel
Temp MAX	-0.5 °C
Temp MIN	-0.5 °C
Wide	520 mm
Height	637 mm
Length	345 mm
Depth	345 mm
Operation via	Knobs/Buttons
Weight	39 kg
Refrigerant medium	Refrigerant (closed system)
Refrigerant closed system	R 290 (propane)
Voltage	230 V
for ambient temperature	10 to 40 °C
Refrigerators and freezers features	suitable for use with deionized water
Ice capacity kg/24h	40 kg/24h
Quantity compressors	1
Capacity reservoir	10 kg
Power	410 W
Devices self cooling	Water cooling
Positioning variants	single
Type refrigerator / freezer	Flake ice maker

