



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Manitowoc flake ice maker, 113 kg/day, air cooling

**Price on
request
plus VAT &
Shipping**

Product Images

Description

Powerful **Manitowoc® brand flake ice maker**
UFP series **with storage**

The units produce different quantities of compact, granular ice flakes particularly quickly and cost-effectively.

Cost-effective and fast production

Flake ice has large surface areas and thus cools particularly quickly. It can be used in a wide range of applications – from supermarkets to food production and laboratories – and placed directly on the products to be cooled. What is always important for all users is the consistently high quality of the compact and granular flakes with very good cooling performance.

Using efficient technologies

Flake ice is produced in compliance with HACCP in a vertical stainless steel evaporator. In this process, a screw conveys the ice in a closed system to the head of the cylinder, where it is pressed and broken into the desired size. In this particularly hygienic and economical process, almost the entire quantity of water used is converted into ice.

Naturally sustainable and highly efficient

All ice makers use the environmentally friendly refrigerant R 290, which protects the environment and reduces energy requirements by up to 30%. The naturally occurring propane gas has no ozone depletion potential and an extremely low global warming potential.

Warranty

Ice maker: 3 years (parts) - 2 years (labor)

Evaporator: 3 years (parts) - 2 years (labor)

Compressor: 5 years (parts) - 2 years (labor)

All devices are suitable for deionized water

Additional Information

| | |
|-------------------------------------|---------------------------------------|
| No. | TW-0010 |
| Manufacturer (Brand) | Manitowoc |
| old neolab article no. | TW-0010 |
| EAN | 4058072973032 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Stainless steel |
| Wide | 500 mm |
| Height | 843 mm |
| Depth | 660 mm |
| Operation via | Knobs/Buttons |
| Refrigerant medium | Refrigerant (closed system) |
| Refrigerant closed system | R 290 (propane) |
| Voltage | 230 V |
| for ambient temperature | 10 to 43° C |
| Refrigerators and freezers features | suitable for use with deionized water |
| Quantity compressors | 1 |
| Ice capacity kg/24h | 113 kg/24h |
| Devices self cooling | Air cooling |
| Interior | Collection container to fold open |
| Capacity reservoir | 20 kg |
| Positioning variants | single |
| Type refrigerator / freezer | Flake ice maker |

