



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Eppendorf CryoCube F740hiw, ultra-low  
temperature freezer, 740L, VIP, WC, DoRi, 5,  
230V/50Hz (EU)**

**€19,080.00  
plus VAT &  
Shipping**

## Product Images

---

## Description

The Eppendorf CryoCube® F740 series is one of the newest additions to the Eppendorf family of ultra-low temperature freezers (also called ultra-low temperature freezers, -80°C). The CryoCube F740 ultra-low temperature freezer series combines the traditional high quality of our previous ultra-low temperature freezers with a higher capacity to safely store more of your samples while optimizing energy consumption. By incorporating the PhysioCare Concept® during development, the CryoCube F740 has been made easier and more comfortable to use, supporting workplace ergonomics and enabling optimal workflow in the lab.

## Additional Information

No.	VB-3987
Manufacturer (Brand)	Eppendorf
EAN	4058072999742
Gross weight	377
HSNumber	84184080
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	328.00000 kg
Country of origin	Great Britain
Refrigerant closed system	R 290 (propane)
Devices self cooling	Water cooling
Capacity for 50 boxes	576
Volume freezer / refrigerator	740 liters

