



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® ultra-low freezer, 100 liters, underbench cabinet

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Ultra-low freezer with low-noise refrigeration unit and 100 liters usable capacity, the ideal unit for storing small quantities of samples or for storage directly at the workstation in the laboratory.

- Two-stage cascade cooling system
- RS485 interface
- Lockable door
- With insertion grid and 4 shelf rails

Technical data:

Temperature range: -80°C to -50°C

Control accuracy: +/- 1K

Housing: stainless steel

Inner dimensions (WxDxH): 450 x 450 x 500 mm

Weight: 120 kg

Ambient temperature: +12°C to +30°C

Order No. 7-4564; Underbench cabinet

Ventilation grille in front and on the right, thus suitable for installation under the table.

Stackable (2 units)

External dimensions (WxDxH): 950 x 730 x 750 mm

Optional on request:

Feed-through 19 mm

4 fixed castors, 2 with locking brake

Additional Information

No.	7-4564
Manufacturer (Brand)	neoLab
EAN	4058072040659
Transport temperature	Room temperature
Color	Grey Blue
Material	Stainless steel
Temp MAX	-50 °C
Temp MIN	-80 °C
Wide	950 mm
Width/inside	450 mm
Height	750 mm
Height inside	500.0 mm
Length	730 mm
Depth	730 mm
Depth inside	450 mm
Number of shelves/ grids	4
Operation via	Membrane keypad
Weight	120 kg
Refrigerant medium	Refrigerant (closed system)
Refrigerant closed system	R 290 (propane)
Length inside	450 mm
Voltage	230 V
Type interface	RS485
for ambient temperature	12 to 30 °C
Refrigerators and freezers features	with potential-free contact for remote alarm with interface with buffer battery for data storage with possibility of execution with data logger with alarm function lockable
Number of implementation options	1
Quantity compressors	1

Number of drawers	0
DM MAX Implementation option	19 mm
Interior	Shelves
Capacity for 50 boxes	8l
Page door suspension	left
Positioning variants	submountable
Type refrigerator / freezer	Ultra Deep Freezer / Deep Freezer
Volume freezer / refrigerator	100 l

