



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Heidolph Hei-TORQUE Core -- EU cable

€1,051.00
plus VAT &
Shipping

Product Images



Description

Heidolph's lightest laboratory stirrer is here - with a torque of 40 Ncm! Thanks to its weight of 2.3 kg, the Hei-TORQUE Core is easier to carry than the other models and to mount on a stand and/or in various closed system installations. However, this does not mean that this stirrer is less powerful. The additional name "Core" emphasizes that this is a model with select functions and important safety features that provide noticeable assistance in any stirring task. Above all, the Hei-TORQUE Core scores with its ease of use: the quick chuck is very easy to open and close without a key. The "Max" button allows you to reach the maximum speed of 2,000 rpm in the shortest possible time. Thanks to the display you keep all application parameters in view. The Hei-TORQUE Core has a quick-action chuck with a diameter of 10.5 mm, allowing it to be combined with a wide range of stirring devices. To achieve the best homogeneity, it is suitable for viscosities up to 10,000 mPas or volumes up to 25 liters (H 2 O).

Additional Information

No.	BH-1276
Manufacturer (Brand)	Heidolph
EAN	4058072774868
Gross weight	3.2
Transport temperature	Room temperature
Net weight	2.80000 kg
Speed MAX	2000 rpm
Country of origin	Germany
Region of origin	Bavaria
Packing width	0.26000 m
Packing height	0.14000 m
Packing depth	0.39000 m
Packaging volume	0.01400000 ccm
Customs tariff number	84798200
Speed MIN	20 rpm
Power	77 W

