



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Heidolph Hei-TORQUE Platinum Package -- EU Cable

**€3,246.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

To provide you with the perfect all-in-one solution for powerful stirring and the easiest operation in the lab, we have added packages to the Hei-TORQUE product group. For ideal use on laboratory benches, each package includes a telescopic stand and a corresponding clamp. Hei-TORQUE Ultimate 400 Telescopic Stand Clamp

## Additional Information

---

No.	BH-1306
Manufacturer (Brand)	Heidolph
EAN	4058072775100
Transport temperature	Room temperature
Speed MAX	2000 rpm
Country of origin	Germany
Region of origin	Bavaria
Customs tariff number	84798200
Speed MIN	10 rpm
Power	150 W

