



Thermo Scientific Variomag Micro Stirrer

Product Images

Description

Variomag Micro Stirrer Work with Thermo Scientific™ Cimarec™ i micro stirrers, IP68 hermetically sealed stirrers, equivalent to protection against continuous submersion Certification/Compliance: Certificate of Analysis (COA), Certificate of Conformity (COC) Description: Micro Stirrer w/o Controller Protection Class: IP68 Type: Stirrer Capacity: 1 L Speed: 1400 RPM Cover Plate Material: Stainless Steel Unit Size: Each

Additional Information

No.	LD-1281
Manufacturer (Brand)	Thermo Scientific
EAN	4066292010604

