



Thermo Scientific Variomag Micro Stirrer

Product Images

Description

Variomag Micro Stirrer Work with Thermo Scientific™ Cimarec™ i micro stirrers, hermetically sealed IP68 rated stirrers, equivalent to protection against continuous submersion Certification/Compliance: CE Description: Micro w / Standard Telemodule; 230 Euro Dimensions (L x W x H): 48 x 48 x 15cm (1.8 x 1.8 x 0.6 in.) Height (anglo): 0.6 in. Height (metric): 15 mm Length (Anglo): 1.8 in. Plug type: EU version Protection class: IP68 Temperature range (metric): -10 ° to 56 °C (operating temp) Type: stirrer Weight (Anglo-American): £ 0.4 Width (Anglo-American): 1.8 in. Capacity: 1 L Shape: Square Speed: 130 - 1000rpm Stirring Positions: 1 Cover Plate Material: Stainless Steel Effective Power: 3w Unit Size: Each

Additional Information

No.	LD-1134
Manufacturer (Brand)	Thermo Scientific
VGKL number	517940107
EAN	4066292009257
Net weight	0.20000 kg
Wide	48 mm
Height	15 mm
Length	48 mm
Voltage	230 V

