



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLab® tumbling roller mixer 5 rollers, speed: 4-30 rpm

**€608.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Synchronously rotating PVC rollers form a uniformly tumbling support surface.

The samples are rotated around their longitudinal axis and simultaneously rocked back and forth along their support surface.

The device is specially designed for uniform and gentle mixing of samples for hematological examinations and other cell suspensions. It prevents sedimentation of samples (dispersions) of all kinds, e.g. paints, varnishes, oils, cosmetics, etc.

Center distance: 40 mm

## Additional Information

No.	ZT-0384
Manufacturer (Brand)	neoLab
old neoLab article no.	3-5075
EAN	4058072444389
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	30 rpm
Color	White
sterile	No
Wide	177 mm
Height	107 mm
Length	445 mm
Depth	445 mm
Weight	3.8 kg
Voltage	230 V
for ambient temperature	5 to 40 °C
Stirrers and rotators features	with speed monitoring with continuous function
Movement	rolling tumbling
Speed MIN	4 rpm
Number of speed levels	stepless
Number of rolls	5
Power	25 W
DM Rollers	32 mm
Maximum relative humidity at 20 °C	80.0 %
MIN tilt angle	3 °

