



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Revolver Rotator, fixed

€92.63
plus VAT &
Shipping

Product Images



Description

The rotator can be used for a variety of mixing tasks in biochemistry, histochemistry, clinical chemistry and molecular biology.

The device has a fixed speed of 18 rpm.

Sample vessels of different volumes can be fixed in any intermediate position.

Both halves of the paddles can be independently positioned between 0-90°.

Fixed speed of 18 rpm.

Easy fixation of the paddles by plug-in connections (without tools).

Paddle for 6 x 50 ml tubes can be used vertically as well as horizontally.

Scope of delivery:

Rotator incl. power supply

2 paddles for 7 x 10/15 ml tubes each

2 paddles for 7 x 5/7 ml blood collection tubes each

2 paddles for each 15 x 1.5/2.0 ml reaction tubes

2 paddles for 21 x 0.5/0.8 ml tubes each

1 paddle for 6 x 50 ml tubes

1 additional rotation axis for paddle 6 x 50 ml

10 spare O-rings

Additional Information

No.	D-6049
Manufacturer (Brand)	neoLabLine
EAN	4058072086916
Transport temperature	Room temperature
Speed MAX	18 rpm
Color	White
Wide	148 mm
Height	195 mm
Length	260 mm
Depth	260 mm
Stirrers and rotators features	with speed monitoring with continuous function
Number of speed levels	1
Attachment is suitable for	Reaction vessels and centrifuge tubes. 7 x 10 to 15 ml (2x) (D-6051) centrifuge tubes. 6 x 50 ml (D-6055) reaction tubes. 15 x 1.5/ 2.0 ml (2x) (D-6053) reaction tubes. 21 x 0.5/ 0.8 ml (2x) (D-6054) Blood collection tubes. 7 x 5.0/ 7.0 ml (2x) (D-6
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes
Operation via	Knobs/Buttons
Movement	rotating vibrating
Speed MIN	18 rpm
Spare part	D-6051 D-6052 D-6053 D-6054 D-6055
Weight	1 kg
Layer tube	variable
Power	4 W
Maximum relative humidity at 20 °C	80.0 %
Model Shaker/ Rotator	neoLab turret rotator (D-6049)
Voltage	230 V
Shaker type	Orbital/circular shaker

for ambient temperature

4 to 60 °C

