



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine rotator adjustable 90° to 27°, 10 rpm fixed, without turntable

€689.90
plus VAT &
Shipping

Product Images



Description

For mixing samples in test tubes or reaction vessels overhead.

Powerful motor with different speeds (fixed or variable)

Power supply cable included.

Power supply: 230 V/50 Hz

Order No. 2-1175

Weight: 2.6 kg

Additional Information

No.	2-1175
Manufacturer (Brand)	neoLabLine
EAN	4058072055721
Transport temperature	Room temperature
Speed MAX	10 rpm
Color	yellow
Wide	200 mm
Height	220 mm
Length	210 mm
Depth	210 mm
Display (reading, setting)	analog
Operation via	Knobs/Buttons
Weight	2.6 kg
Voltage	230 V
Stirrers and rotators features	with continuous function turntables are available as accessories
Movement	rotating
Speed MIN	10 rpm
Number of speed levels	1
Power	6 W
Model Shaker/ Rotator	neoLabLine Rotator 10 rpm fix (2-1175)
MAX tilt angle	90 °
MIN tilt angle	27 °
Shaker type	Rotator/ Overhead Shaker

