



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Rotator, 5-50 rpm, with digital display, variable, without turntable

€1,201.00
plus VAT &
Shipping

Product Images



Description

For mixing samples in test tubes or reaction vessels overhead.

Powerful motor with different speeds (fixed or variable)

Power supply cable included.

Power supply: 230 V/50 Hz

Order No. 2-1185

Digital display for speed and time (0 - 99 h / 1 - 59 min), settings can be saved.

Weight: 3.8 kg

Additional Information

No.	2-1185
Manufacturer (Brand)	neoLabLine
EAN	4058072055813
Transport temperature	Room temperature
Speed MAX	50 rpm
Color	yellow
Wide	200 mm
Height	250 mm
Length	200 mm
Depth	200 mm
Display (reading, setting)	digital
Operation via	Membrane keypad
Weight	3.5 kg
Voltage	230 V
Stirrers and rotators features	with infinitely variable speed control with timer with continuous function rotary plates are available as accessories
Number of speed levels	stepless
Movement	rotating
Speed MIN	5 rpm
Power	50 W
Model Shaker/ Rotator	neoLabLine Rotator 5 to 50 rpm digital (2-1185)
MAX tilt angle	90 °
MIN tilt angle	27 °
Shaker type	Rotator/ Overhead Shaker

