




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine Rotator with vortexer, without rack

€680.00
plus VAT &
Shipping

Product Images



Description

Variable speed rotator combined with the advantages of a vortexer.

From gentle hybridization at low speed to fast vortexing, everything is possible thanks to the different racks and operating variants.

- 34 different programs
- 6 variable vortexer programs
- 3 additional individually storable programs
- Maintenance free motor
- Combi Rack can be equipped with 2 x 50 ml, 4 x 15 ml, 8 x 5 ml and 8 x 1.5 ml

Technical data

Rotation: 1 to 99 RPM (increments of 1)

Programs: 34

Ambient temperature: 2 to 50°C

Rel. humidity: max. 80 %

Dimensions (WxDxH): 344x125x168 mm

Weight: 1.6 kg

Mains connection: 12 V, incl. AC adapter, 1500 mA

Accessories: Various rack attachments (ArtNo. 7-0043 and 7-0044, 7-0046 to 7-0050, 7-0053)

Additional Information

No.	7-0045
Manufacturer (Brand)	neoLabLine
EAN	4058072026752
Transport temperature	Room temperature
autoclavable at 121°C	No
Speed MAX	99 rpm
Color	yellow
sterile	No
Wide	125 mm
Height	168 mm
Length	344 mm
Depth	344 mm
Stirrers and rotators features	programmable with infinitely variable speed control with continuous function
Number of speed levels	stepless
Attachment is suitable for	Blood Collection Tubes. 28 x 5 ml (7-0048) reaction tubes. 32 x 0.5 ml (7-0044) reaction tubes. 32 x 1.5/ 2.0 ml (7-0049) centrifuge tubes. 10 x 50 ml (7-0046) Reaction vessels and centrifuge tubes. 8 x 1.5/2.0 ml. 4 x 2.2 ml. 4 x 5/15 ml. 2 x 50 ml
Attachment is suitable for type vessel	Reaction vessels cryogenic tubes centrifuge tubes test tubes/tubes
Operation via	Membrane keypad
Movement	teetering
Speed MIN	1 rpm
Weight	1.6 kg
Power	12 W
Maximum relative humidity at 20 °C	80.0 %
Model Shaker/ Rotator	neoLab rotator with vortexer (7-0045)
Voltage	12 V
Shaker type	Rotator/ Overhead shaker Vibrating shaker/ Vortexer Rocking shaker

for ambient temperature

2 to 50 °C

