



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® USB Microscope, 50–500x

**Price on
request
plus VAT &
Shipping**

Product Images



Description

These compact USB mini microscopes are the simple solution for quick and easy optical analysis of small parts and other microworlds. Whether in the lab, in the workshop or in the open field, the lightweight portable devices allow an examination in any place. 8 LEDs provide a bright and optimal display of the things to be examined. The devices are compatible with Android, PC, Mac and Linux devices.

7-3160: Resolution 50-500x

Additional Information

No.	7-3160
Manufacturer (Brand)	neoLab
EAN	4058072122188
Transport temperature	Room temperature
Model microscope	USB Microscope for Windows / Mac-OS (x) / Linux
Type connection cable	USB
Type illuminant	LED
Type microscope	USB Microscope
Type tube camera	monocular
external interface camera	USB 2.0
with dust cover	false

