



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN Metallurgical Microscope (Inverted) Monocular Achromat 10/20/40: WF10x18: 3W LED (IL)

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

- The KERN OLE-1 is a very easy to use, stable inverted metallurgical microscope
- It is suitable for all common routine applications in training facilities or workshops and provides excellent images
- Available in monocular or binocular versions with different eyepieces
- As standard, these microscopes are equipped with planachromatic objectives
- A revolver for up to 4 objectives and a large specimen stage are also included as standard
- The following optional accessories are available: LWD objectives for a long working distance as well as different stage inserts etc.
- A powerful 3W LED illumination is available for illumination
- A central feature of this variable and at the same time robust microscope series is the stable and precisely adjustable mechanics. This is underlined by the functional and ergonomic design

## Additional Information

---

No.	KP-5042
Manufacturer (Brand)	KERN
EAN	4058072145965
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Upright microscopes
Manufacturer	CORE
Inclination tube	36 °
Type lighting	External light source
Type tube camera	monocular
with dust cover	false

