



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Zeiss Axiolab 5 binocular microscope for transmitted-light brightfield and phase contrast with phototube

€7,783.00
plus VAT &
Shipping

Product Images



Description

Preconfigured microscopes from Zeiss, specially selected for you by neolab.

This is an example configuration which we will be happy to adapt to your personal needs and requirements. Please contact our product specialists Dr. Hahn or Jochen Weber directly.

The **Zeiss Axiolab 5** is designed for routine microscopy performed daily in your laboratory. The compact and ergonomic design saves space and allows easy handling. Axiolab 5 is a true team player. Combine it with Axiocam 208 color and benefit from the intelligent microscopy concept: you will experience a completely new form of digital documentation. Simply focus your specimen and press a single button for razor-sharp images in true color. The digital image looks as if you are looking at it through the eyepieces, with all the details and subtle color differences clearly visible.

In addition, Axiolab 5 automatically adds the correct scaling information to your images. You get all this in standalone mode, without needing a PC or additional software. With Axiolab 5 you can save time, money and valuable lab space. Digital documentation has never been easier.

Application area: cytology, microbiology, gynecology, cytology

Typical samples : Body fluids, gynecological smears, bacteria, viruses, fungi, parasites, semen

Please note that full functionality of a pre-configured instrument cannot be guaranteed if an item is removed or added.

Additional Information

No.	CZ-0520
Manufacturer (Brand)	Zeiss
old neoLab article no.	7-6631
EAN	4066292499959
Transport temperature	Room temperature

