



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN stereo microscope OSF, 37653

€410.00
plus VAT &
Shipping

Product Images



Description

- The KERN OSF-4G is specially designed for schools and workshops due to its integrated handle as well as its stable mechanical stand
- The LED incident light and transmitted light illumination included as standard ensures continuously dimmable and optimal exposure of your specimen
- In addition to the very good optical properties, it offers the highest comfort in this class due to its ergonomically designed work surface
- An interchangeable objective with predefined magnifications is available for fast and efficient work
- The eyepieces are fixed in the tube to protect them from loss or damage
- The ergonomic shape and the stable and extremely precisely adjustable mechanics offer high functionality and enable fast and efficient work with only a few hand movements
- A wide range of eyepieces and various additional external illumination units are available as accessories
- A dust cover, eyecups, as well as a multilingual operating manual are included in the scope of delivery
- KERN OSF 439: Also available in a set of 10

Field of application: Education, in vitro fertilization, detection of parasites, zoology and botany, tissue preparation, dissection, quality control

Applications/specimens: preparations with focus on spatial impression, zoom with variable magnification (depth, thickness) e.g. insects, seeds, circuit boards, components

Technical data:

- Optical system: Greenough optics
- Illumination independently dimmable
- Tube tilted 45°
- Interpupillary distance: 55.0-75.0 mm
- Dioptric compensation on one side
- Overall dimensions W×D×H 230×180×275 mm
- Net weight approx. 2.5 kg

Additional Information

No.	KP-5131
Manufacturer (Brand)	KERN
EAN	4058072146801
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Manufacturer	CORE
Type tube camera	binocular
with dust cover	false

