



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereo microscope OSF, 43525

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

- The KERN is a very robust, stable and easy to use stereo microscope
- Due to its design and the flexible column stand, it is ideally suited for rough working environments, workshops, schools and training centers
- In addition to its very good optical properties, it offers the highest level of comfort in this class due to its large work surface
- The powerful and separately dimmable incident and transmitted light LED illumination unit, which can be selected depending on the model, always reliably provides good illumination
- The variety of objective combinations, the extra-long working distance and the depth of field leave nothing to be desired
- The large selection of additional eyepieces, also for eyeglass wearers, as well as other accessories, rounds out this package
- A central feature of this variable and at the same time robust microscope series is the stable and precisely adjustable mechanics. This is underlined by the functional and ergonomic design

## Additional Information

No.	KP-5138
Manufacturer (Brand)	KERN
EAN	4058072146870
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type microscope	Stereomicroscopes
Manufacturer	CORE
Type tube camera	binocular
with dust cover	false

