



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereomicroscope set OZL, 0.7 x - 4 5

**€1,090.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- The KERN OZL-96 belongs to the stereo zoom microscopes that convince with their quality, easy handling, flexibility as well as their stability and favorable price
- The 4.5W LED ring illumination included as standard ensures optimal illumination of your specimen
- In addition to the good optical properties, these models offer the highest comfort in this class due to their large work surface - optimal for training facilities, workshops as well as assembly and repair workplaces, e.g. in the electronics industry
- The zoom lens allows you to magnify continuously from 7×-45×
- The KERN OZL-9 series is available as binocular or trinocular version
- The universal telescopic arm stand offers you a high degree of flexibility. Especially suitable for large specimens and free working under the microscope
- Pre-defined sets consisting of a stereo microscope head, a universal stand and a ring illuminator from our range of products
- This saves you the time-consuming configuration and the agony of choosing the combination of the various components. So you get an inexpensive and highly flexible solution for your microscope workstation
- Simple - practical - inexpensive

## Additional Information

---

No.	KP-5266
Manufacturer (Brand)	KERN
EAN	4058072148027
Transport temperature	Room temperature
Energy supply	Power grid
Manufacturer	CORE
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Type tube camera	trinocular
with dust cover	false

