



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN stereomicroscope set OZM, 0.7 x - 4 5

€1,750.00
plus VAT &
Shipping

Product Images



Description

- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 provides true color and depth sharp images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The extendable telescopic arm stand offers you the possibility to adjust the height individually to your specimen. An additional extendable column in the horizontal axis of the telescopic arm guarantees you an additional reach

Additional Information

No.	KP-5279
Manufacturer (Brand)	KERN
EAN	4058072148119
Transport temperature	Room temperature
Energy supply	Power grid
Manufacturer	CORE
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Type tube camera	binocular
with dust cover	false

