



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereomicroscope set OZM, 0.7 x - 4 5

**€2,090.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 provides true color and depth sharp images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The extendable telescopic arm stand offers you the possibility to adjust the height individually to your specimen. An additional extendable column in the horizontal axis of the telescopic arm guarantees you an additional reach

## Additional Information

---

No.	KP-5280
Manufacturer (Brand)	KERN
EAN	4058072148126
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Manufacturer	CORE
Type tube camera	trinocular
with dust cover	false

