



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereomicroscope set OZM, 0.7 x - 4 5

**€2,080.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 delivers true color and depth of field images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The ball bearing double arm universal stand gives you smooth, precise and easy adjustment in the horizontal axis

## Additional Information

No.	KP-5286
Manufacturer (Brand)	KERN
EAN	4058072148171
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Manufacturer	CORE
Type tube camera	binocular
with dust cover	false

