



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN stereomicroscope set OZM, 0.7 x - 4 5

€2,420.00
plus VAT &
Shipping

Product Images



Description

- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 delivers true-color and depth-sharp images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The ball bearing double arm universal stand gives you smooth, precise and easy adjustment in the horizontal axis

Additional Information

| | |
|-----------------------|-------------------|
| No. | KP-5288 |
| Manufacturer (Brand) | KERN |
| EAN | 4058072148195 |
| Transport temperature | Room temperature |
| Energy supply | Power grid |
| Physical principle | Transmitted light |
| Type illuminant | LED |
| Type microscope | Stereomicroscopes |
| Manufacturer | CORE |
| Type tube camera | trinocular |
| with dust cover | false |

