



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN stereomicroscope set OZM, 0.7 x - 4 5

€2,090.00
plus VAT &
Shipping

Product Images



Description

- Already predefined sets consisting of a stereo microscope head, a universal stand and a ring illuminator from our range
- This saves you the time-consuming configuration and the agony of choosing the combination of the various components. So you get an inexpensive and highly flexible solution for your microscope workstation
- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 delivers true color and depth sharp images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The individually adjustable articulated arm universal stand with clamp guarantees you the highest flexibility and reach. Due to the individual height adjustment and extendable articulated arms, large specimens can be examined
- The clamp allows the universal stand to be attached to any work table
- Simple - practical - inexpensive

Additional Information

No.	KP-5293
Manufacturer (Brand)	KERN
EAN	4058072148249
Transport temperature	Room temperature
Energy supply	Power grid
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Manufacturer	CORE
Type tube camera	trinocular
with dust cover	false

