



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN stereomicroscope set OZM, 0.7 x - 4 5

**€2,210.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

- The KERN OZM-9 includes outstanding stereo zoom microscopes with above average optical performance
- The 4.5W LED ring illuminator included as standard ensures optimal illumination of your specimen
- In addition to the long working distance, an extra large field of view and its brilliant resolution, the KERN OZM-9 delivers true color and depth sharp images with high contrast
- The zoom objective allows you to magnify continuously from 7×-45×
- You can choose between a binocular model and a trinocular model for connecting a camera for documentation and quality reports
- The gas spring articulated arm with clamp allows you to make quick and easy adjustments. The set position is held by the built-in gas spring

## Additional Information

---

No.	KP-5296
Manufacturer (Brand)	KERN
EAN	4058072148270
Transport temperature	Room temperature
Energy supply	Power grid
Manufacturer	CORE
Physical principle	Transmitted light
Type illuminant	LED
Type microscope	Stereomicroscopes
Type tube camera	trinocular
with dust cover	false

