



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Zeiss microscope body Stemi 305 trino for retrofitting a camera

**€996.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Zeiss microscope body Stemi 305 is a stereo microscope based on the Greenough principle and delivers images of the highest brilliance in unique quality. With a working area of 195 x 160 mm with a working distance of 110 mm, a wide range of applications is made possible. Supplements

- Stereomicroscope with 45° binocular tube
- Manual, ambidextrous zoom 5:1 (0.8x to 4.0x)
- Fixed zoom detent positions 0.8x-1x-2x-3x-4x
- Vertically integrated LED illumination (IVI) for tripods K EDU/LAB/MAT or controller K LED
- Eyepieces 10x/23 Br. foc.
- Spiral cable RJ12, tripod K
- Built-in power supply 12V DC 24W/100...240V AC/50...60Hz
- Incl. power cable EURO C7 and dust cover

Order No. 7-1006: Microscope body Stemi 305 trino for retrofitting a camera

## Additional Information

---

No.	CZ-0477
Manufacturer (Brand)	Zeiss
old neoLab article no.	7-1006
EAN	4058072226640
Transport temperature	Room temperature
Number of eyepieces	2
Physical principle	Polarization and attachment optics. optional. 1 connector (M52)
Type illuminant	LED
Type microscope	Stereomicroscopes
Type eyepiece	10x/23 mm
Type lighting	Incident light
Type tube camera	trinocular
with dust cover	false

