



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND BioChem-VacuuCenter BVC control G, with 2l glass bottle, 230 V / 50-60 Hz, CEE Net

€1,945.00
plus VAT &
Shipping

Product Images



Description

Liquid suction system BVC control G --- --- Max. Suction capacity 50/60 Hz (m³/h): 0.7/0.85 --- Ultimate vacuum (abs.) (mbar): 150 --- Vacuum controller: touch panel --- Connection suction side: suction handle --- Connection discharge side: hose shaft DN 8 mm / sound insulation. --- Dimensions (L x W x H), approx. (mm): 408x194x430 --- Weight, approx. (kg): 7.7 --- Rated voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Rated motor power (kW): 0.04 --- Scope of delivery: Liquid suction system complete with powerful chemistry diaphragm pump, collecting bottle, suction power control adjustable via touch panel, suction handle VHCpro and protective filter, ready for use with instructions. Old item number: 727300

Additional Information

No.	VZ-1466
Manufacturer (Brand)	VACUUBRAND
VGKL number	101877273
EAN	4066292394919
Gross weight	11.5
Transport temperature	Room temperature
Net weight	8.30000 kg
Country of origin	Germany
Customs tariff number	84141089

