



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



VACUUBRAND BioChem-VacuuCenter BVC professional, with 4l PP bottle, 230 V / 50-60 Hz, CEE N

€2,430.00
plus VAT &
Shipping

Product Images



Description

Liquid suction system BVC professional --- --- Max. Suction capacity 50/60 Hz (m³/h): 0.7/0.85 --- Ultimate vacuum (abs.) (mbar): 150 --- Vacuum controller: touch panel --- Connection suction side: suction handle --- Connection pressure side: hose shaft DN 8 mm / sound insulation. --- Dimensions (L x W x H), approx. (mm): 408x194x500 --- Weight, approx. (kg): 7.3 --- Rated voltage / mains frequency: 230 V/50-60 Hz --- Plug: CEE --- Motor rated power (kW): 0.04 --- Scope of delivery: Liquid suction system complete with powerful chemistry diaphragm pump, collection bottle, suction power control adjustable via touch panel, non-contact level sensor, self-closing quick couplings (for 4l PP bottle), suction handle VHCpro and protective filter, ready for use with instructions. Old item number: 727400

Additional Information

No.	VZ-1470
Manufacturer (Brand)	VACUUBRAND
VGKL number	101877274
EAN	4066292394957
Gross weight	11.5
Transport temperature	Room temperature
Net weight	8.30000 kg
Country of origin	Germany
Customs tariff number	84141089

