



Asecos venting module

Product Images



Description

Ventilation top W x D x H (mm): 200 x 400 x 200 Information/Group: with exhaust air monitoring Material/Color/Description: sheet steel powder-coated RAL 7035 Model: suitable for cabinets S-LINE / Q-LINE / G-LINE* old item number: HF.EA.11442. Delivery type: incl. free delivery Germany (mainland)

Additional Information

No.	AS-0197
Manufacturer (Brand)	Asecos
EAN	4058072479695
Net weight	10.10000 kg
Country of origin	Deutschland
Region of origin	Hessen
Customs tariff number	84145100

