



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## VACUUBRAND suction system with 4 ltr. PP bottle

**€895.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The suction system is designed for connection to an existing external vacuum source such as a diaphragm pump or vacuum network and does not require an electrical supply.

The chemical-resistant, mechanical vacuum regulator ensures an optimum working vacuum. The regulator opens only when needed, limiting the gas load for the vacuum source.

Evaporation of the collected liquid is minimized. The guide tube integrated in the bottle head serves to reduce foam and aerosol formation and ensures long service life of the hydrophobic protective filter.

### Performance features

Mechanical vacuum regulator to minimize gas load for the vacuum network and for optimum working vacuum. High stability, no electrical connection required for flexible use. Optimally equipped autoclavable collection bottle to reduce aerosol and foam formation, autoclavable quick disconnect couplings available as accessories. Suction handle with continuous suction hose - mechanics not in contact with media - no contamination. Hydrophobic protective filter 0.2 µm - high safety against contamination of the vacuum network, autoclavable as a unit with the collection bottle.

Weight: approx. 2.8 kg

Accessories: Spare collection bottle (Order No. 7-8156 and 7-8157)

Quick coupling complete (Order No. 7-8158)

## Additional Information

No.	VZ-1459
Manufacturer (Brand)	VACUUBRAND
old neoLab article no.	7-8154
EAN	4058072426743
Transport temperature	Room temperature
Wide	200 mm
Height	490 mm
Length	250 mm
Depth	250 mm
Weight	2.8 kg
Pumps and suction systems Features	with safety filter with collection bottle
Functionality	mechanical
for DM hose inside 1	6.0 to 10.0 mm
Spare part	VZ-0116
Material collection bottle	Polypropylene (PP)
Type suction device/system	without pump
Type connection collection bottle	Hose olive
Volume collection bottle	4 l

